# PREZA-PAK II

## ARTERIAL BLOOD SAMPLER

### ■ INSTRUCTIONS FOR USE

1. Remove rubber cap and attach blood collecting needle. (\*1)
2. Move plunger and set the gasket tip to the graduation mark necessary to obtain the blood volume required. (Fig. 1)
3. Disinfect the site and perform arterial puncture with the sampler at 45°—90° angle from horizontal. (\*2) (Fig. 2)
4. After completing blood collection, immediately place into the seal protector.
5. Thoroughly mix heparinized sample by rolling between palms about 50 times.
6. Transport to blood analysis area for processing.
  - \*1. Use 23G or 22G disposable needle for blood collecting needle. When less than 1 ml blood volume is required, use 23G or 25G disposable needle.
  - \*2. Take care that blood does not touch the filter before blood collection is completed, otherwise air in the sampler will not come out.

#### Two Alternative Method For Aspiration

1. When normal collecting procedure above does not allow sufficient blood in flow, close the opening at the end of plunger with a finger and pull the plunger slowly. (Fig. 3)

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2. When the opening at the end of plunger is closed with the rubber cap, pull the plunger and collect blood through aspirating straightly in the same way as blood collection with a syringe. (Fig. 4)

#### How To Remove Remaining Air In Sampler:

Caution: During the following procedure, PREZA-PAK should always be kept vertically, after collecting blood.

1. Put the protector immediately after blood collection.
2. Make the opening at the end of plunger loose and push the plunger down slowly.

#### Note

1. For in-vitro diagnostic use (blood gas analysis).
2. If the sample is not analyzed immediately, keep in refrigerator.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

◆ Trademark

©, TERUMO and PREZA-PAK are Registered in U.S. Patent and Trademark Office by TERUMO CORPORATION.

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